

DMC-EZ 2800

DMC-EZ

DICOM Media Creators

Product Overview

DMC-EZ Medical Software

- Intuitive Browser-based User Interface
- Creates IHE DICOM-compliant CDs and DVDs
- On-demand CD/DVD Creation for DICOM Studies
- Automatic Disc Creation using Powerful Rules Based Scheduling
- Multiple AE Title Capability
- Includes Multiple Viewers for Disc Distribution
- Include Diagnostic Reports and Non-DICOM Data on CD
- Ability to View Complete Report on Disc with Images and Hyperlinks
- Flexible Disc Labeling Options
- Study Disc Anonymization
- View Thumbnail Images prior to Disc Creation
- Built in Web Server and Web Viewer for Alternative Distribution of Images and Reports
- Image Distribution via USB Flash for Easy Use by Referring Physicians
- Study Disc Encryption
- Rules Based DICOM Routing
- All CD/DVD Disc Activity Recorded in Audit Logs (HIPAA)
- Back-Up Configuration and Database to USB or CD/DVD Media Which Allows Easy Restore
- Import Data from CDs and DVDs Your Facility Receives

Rimage 3400 Hardware

- Fully Automatic Disc Publishing System
- Two CD/DVD Recorders
- Everest 400 Thermal Re-Transfer Printer
- Scratch proof, water proof, no fade print technology
- Three 50 Disc Input/ Output bin / 5 Disc External output

DMC-EZ PC

- High performance CPU
- Monitor, Keyboard and Mouse



R I M A G E®

- **Small to Moderate Volume System**
- **Small to Medium Imaging Centers and Radiology Departments**
- **Single Modalities or Modality Clusters**
- **50-150 Discs Per Day**



Sales Support: 800.944.6281
Outside US 408.350.3450
fax: 408.350.3100

309 Laurelwood Rd. #20
Santa Clara, CA 95054
www.sstgroup-inc.com
sales@sstgroup-inc.com

DMC-EZ 2800

Features and Benefits

- **Intuitive Browser-based User Interface**

Easy-to-use, intuitive user interface, local language support, designed in cooperation with leading practitioners.

- **DICOM Compliant**

The DMC-EZ works in any standard DICOM environment, allowing hospitals and clinics to simply burn and label a customized set of DICOM studies onto CD/DVD. The CD/DVD study data can be read by any DICOM compatible viewer or imported into any DICOM-compatible system/archive. The CD/DVD media created by the DMC-EZ are fully IHE compliant.

- **Automated and/or Manual Operating Modes**

The DMC-EZ has the ability to automatically create CDs and DVDs using rule-based schedules or manually create media using our web interface. The rule-based schedules along with the image cache on the DMC-EZ server allow the facility to optimize their robotic hardware to handle large and small jobs. The DMC-EZ software can be configured to run bulk media writes at night as well as intelligent spanning of data across discs to minimize frustrating "disc swapping."

- **Multiple AE titles and Scheduling**

DMC-EZ systems have the capability to maintain multiple AE titles, and with our rule-based scheduling, can be configured to automate most of the disc creation workflow. An example would be a rule that automatically creates a CD for each CT exam sent to the DMC-EZ. (AE title and corresponding schedule would be "CT Patient CD") In addition, we could implement a date/time schedule which would archive all of the CT Studies to DVD at the end of the day. Media creation rules can be easily configured using AE Title, Sending Device or DICOM Tag.

- **DICOM Viewer Included**

Discs created with DMC-EZ systems include a series of industry standard DICOM viewers (e.g. EZDICOM, OsiriX - Macintosh viewer) for clients without their own DICOM viewing software. This makes all your DICOM study discs "ready to read".

Many OEM viewers can be added, please contact sales if you have questions.

- **Include Diagnostic Reports and Non-DICOM Data on CD**

DMC-EZ can fetch reports from PACS, PACS Broker and RIS utilizing SR or HL7. Software can also attach Non DICOM Data to studies (i.e. pdf, word files, etc..)

- **Ability to View Complete Report on Disc with Key Images and Hyperlinks**

DMC-EZ also has the ability to integrate the Diagnostic Report with Key Images (Hyperlinks), leading referring physicians straight to the most relevant images.

- **Labeling Capabilities**

The DMC-EZ medical recording software allows for simple generation of unique and customizable labels including logos, patient demographics, volume and study information as well as customized user comments.

DMC-EZ 2800

Features and Benefits

- **Study Disc Anonymization**

All DMC-EZ models allow anonymization of discs to ensure patient confidentiality of images used in instructional and conference presentations.

- **View Image Thumbnails Prior to Disc Creation**

DMC-EZ gives the user the ability to view series and image details prior to disc creation. Features help eliminate costly errors and saves valuable time at each site.

- **Image Distribution via Web Viewer for Easy Use by Referring Physicians**

The Easy-to-use web viewer is designed to allow practitioner to view studies with reports in seconds. Full featured reference viewer provides a complete set of Radiology measurement and annotation tools (Distance, Angle, ROI, Reference Lines, etc.). Multi-modality support: single frame (CT, MRI, CR) and multi-frame (U/S, XA) capability. User can easily print selected images to Windows printer for reference.

- **Image Distribution via USB Flash for Easy Use by Referring Physicians**

DMC-EZ systems can create zip files of patient studies and record to USB Flash devices. Zip files can also be sent via email to referring physicians.

- **Encryption**

DMC-EZ software provides the option to encrypt data on the disc. This has become increasingly important in some countries and states with stringent privacy concerns and need for regulatory compliance.

- **DICOM Routing**

DMC-EZ includes rules based DICOM routing as a standard feature. Automated routing allows the facility to send exams to a specific destination based on unique DICOM tags.

- **HIPAA Readiness**

All DMC-EZ activity is recorded in audit logs and archive logs. Logs can be exported via email or network transfer on an automatic basis each day. These logs enable the facility to track CD/DVD activity for HIPAA purposes and also can provide a shelf management system for media stored in house.

- **Automated Back-Up and Restore**

DMC-EZ provides a simple utility which can be set to automatically back up all configuration, label, database and network settings to USB Media, CD/DVD or Network. The utility provides enhanced serviceability as the backup provides an easy way to restore the system if hardware is replaced.

- **Media Import from DICOM Discs**

DMC-EZ systems include an import utility "EzDicom Sender" which makes it easy to import data from DICOM Part 10 compatible CD/DVDs received from outside facilities. The user can simply import selected patients and/or studies and send to PACS Server or other DICOM destination.

DMC-EZ 2800

Hardware Overview

- **Rimage 3400 Professional System**

High performance dual CD/DVD recorder thermal disc publishing system from Rimage. Robotic print and CD/DVD recording speeds have been optimized for fastest possible throughput.

Thermal Re-Transfer Printing

- Everest 400 color thermal re-transfer printer with high resolution print capability. Scratch proof, water proof and no fade print technology.

Fully Automated Disc Creation

- The 3400 disc publishing system offers higher reliability and performance than typical desktop products. This system includes two CD/DVD recorders, three 50 disc input/output bins and 5 disc side output.

Specifications

Rimage 3400

3400 is an economical, dual recorder solution designed for medium volume applications. A small, compact design accommodates any desk or countertop. Full color or monochrome discs are printed with a consistent cost and high quality using the Everest 400 printer. Features dual CD/DVD front swappable recorders, Three 50 disc input/output bins, 5 disc external output.

Power: 100VAC - 240VAC $\pm 10\%$, 60/50Hz

Operating Temp: Temperature: 60° - 86°F (16°C - 30°C)

Network: 10/100/1000 Ethernet interface

Physical Specifications -

14.8" W x 19" H x 19.5" D (375.92mm W x 482.6mm H x 495.3mm D)

64 lbs. (29Kg)

DMC-EZ PC

Power: 100 - 240VAC $\pm 10\%$, 50/60Hz

Operating Temp: Temperature: 60° - 86°F (16° - 30°C)

Network: 10/100/1000 Ethernet interface

Physical Specifications -

7.36" W x 16.10" H x 17.52" D (186mm W x 409mm H x 445mm D)

41 lbs. (18.5Kg)

- **Small to Moderate Volume System**
- **Small to Medium Imaging Centers and Radiology Departments**
- **Single Modalities or Modality Clusters**
- **50-150 Discs Per Day**