



MEDICAL
RECORDING • SOLUTIONS

DICOM MOD Importer

Simplified importing of DICOM MOD stored images to your PACS

The DICOM MOD Importer allows you to easily import existing DICOM MOD (Magneto Optical) images (Windows or MAC based) and archives to your PACS. With our exclusive USB MOD drive you can use virtually any PC on your network with a USB 2.0 port.

The DICOMReader™ software will allow images to be automatically sent directly to your PACS server. You can establish rules to check patient demographics against your worklist or your PACS server and map changes before they are sent to your PACS server.

Getting your patient data off the shelf and on to your PACS can be as simple as a single click with the SST DICOM MOD Importer solution powered by DICOMReader software.



SST Group Inc.
309 Laurelwood Rd.
Suite 20
Santa Clara, CA 95054

408-350-3450
800-944-6281
408-350-3100 (fax)
www.sstgroup-inc.com

Features (Highlights)	DICOMReader	DICOMReader PRO
Read in most DICOM files from MOD	X	X
Full Media scan	X	X
Easy user interface	X	X
Tree Based study selection process		X
Decodes compressed images	X	X
Multiple DICOM Destinations	X	X
Auto send to Default Destination	X	X
Auto send to Destination (rules based)	X	X
Auto Creation of accession number	X	X
DICOM-WORKLIST client for integrated workflow		X
Demographic verification with PACS		X
Enhanced change of DICOM tags		X
Compliant to Medical Reg. 93/42/EEC	X	X

Models

Part#

DICOM MOD Importer DIMMODSTD
(Includes T5200u, USB cable, power cord, DICOMReader software)

DICOM MOD Importer PRO DIMMODPRO
(Includes T5200u, USB cable, power cord, DICOMReader PRO software)




Other MOD Products

SST MOD Drives - T5200u/T5200i/T5200e/T9100u/T9100i/T9100e

MOD Media - 600MB/650MB/1.2GB/1.3GB/2.3GB/2.6GB/4.1GB/4.8GB/5.2GB/9.1GB

Call SST Group for all of your Magneto Optical requirements



408-350-3450
800-944-6281
www.sstgroup-inc.com