



# DICOM Reader™



Simplified importing of DICOM CD/DVD stored images to your PACS



The DICOM Reader software allows you to easily import existing DICOM CD/DVD images and archives to your PACS.

You can use virtually any PC on your network.

The DICOMReader™ software will allow images to be automatically sent directly to your PACS server. You can establish rules to check patient demographics against your worklist or your PACS server and map changes before they are sent to your PACS server.

Getting your patient data off the shelf and on to your PACS can be as simple as a single click with the DICOMReader software.

[www.sstgroup-inc.com](http://www.sstgroup-inc.com)  
800-944-6281

Features (Highlights)	DICOMReader	DICOMReader PRO
Read in most DICOM files from CD/DVD	X	X
Full Media scan	X	X
Easy user interface	X	X
Tree Based study selection process		X
Decodes compressed images	X	X
Multiple DICOM Destinations	X	X
Auto send to Default Destination	X	X
Auto send to Destination (rules based)	X	X
Auto Creation of accession number	X	X
DICOM-WORKLIST client for integrated workflow		X
Demographic verification with PACS		X
Enhanced change of DICOM tags		X
Compliant to Medical Reg. 93/42/EEC	X	X


Part#

DICOM Reader

DICOM Reader PRO



MEDICAL  
RECORDING • SOLUTIONS



SST Group Inc.  
309 Laurelwood Rd. #20  
Santa Clara, CA 95054  
800-944-6281  
408-350-3100 (fax)  
www.sstgroup-inc.com