



309 Laurelwood Rd. Suite 20
Santa Clara, CA 95054
sales@sstgroup-inc.com

Sales Support: 800.944.6281
Outside US: 408.350.3450
Fax: 408.350.3100

DMC - EZ Hybrid e400

Best of Both Worlds- CD/DVD and Cloud

The DMC-EZ™ Hybrid e400 is the next generation in DICOM image delivery systems. It is an innovative, multi-function appliance designed to give healthcare facilities the ability to deliver DICOM images in a multitude of formats. From hard copy DICOM CD/DVDs to secure electronic delivery of images with a viewer, the DMC-EZ Hybrid e400 addresses all of your DICOM image delivery needs.



DMC-EZ Hybrid Medical Software

- Intuitive browser-based user interface.
- Compatible with large format studies (TOMO, Cardio and Ultrasound)
- Creates IHE DICOM-compliant CDs and DVDs.
- On-demand CD/DVD creation for DICOM studies (Query/Retrieve).
- On-demand cloud delivery of DICOM studies (HITECH compliant)
- Automatic disc creation using powerful rule-based scheduling utilizing multiple AE Titles.
- Includes multiple viewers for disc distribution. Windows viewer supports XP through Windows 10.
- Includes diagnostic reports on CD via HL7 feed, SR, Mitra Broker. We also integrate with major RIS/PACS OEM to obtain diagnostics reports. Contact sales for the options.
- View complete report on disc with images and hyper-links.
- Flexible disc labeling options.
- Study disc anonymization.
- View thumbnail images prior to disc creation.
- Built-in web server and web viewer for alternative distribution of images and reports.
- Optional REST API allows burner to be integrated into your workflow.
- Image distribution via USB flash for easy use by referring physicians.
- Automated or Manual study disc encryption (HIPAA/HITECH 256 AES/FIPS)
- Manual and Automated Rules-based DICOM routing.
- Robot supplies status monitor
- All activity recorded in audit logs (HIPAA Compliance)
- Back-up configuration and HIPAA logs to USB, CD/DVD or network location to allow easy restore.
- Import data from CDs and DVDs your facility receives (EzDicom Sender).
- Multiple language support - interface and viewer (English, French, Italian, Portuguese, Spanish).
- User login supports active directory.
- **Patented technology.**

Epson PP100ii Discproducer

High performance dual CD/DVD recorder inkjet desktop disc publishing system with vibrant six-color inkjet printer and precision AcuGrip™ technology , automated robotic print and CD/DVD recording.

DMC-EZ Hybrid Server and Software

High performance CPU running Windows® 7 Pro or Windows® 10 Pro contains the burning and printing software which controls Epson Discproducer , DMC-EZ hybrid application and image cache. Connects to DICOM network and/or DICOM workstation(s).

Virtual Media Distribution

Simple, fast and secure method to get DICOM images to your patients, referring physicians and other facilities. There are no per study fee's, no special sign-up or accounts to setup. By selecting the study(s), entering an email address and pressing "OK" your facility can send and store up to 400 (1 GB) studies in the cloud at any given time. Each with it's own expiration date. No management of the cloud based data is required. Recipient receives an email with a download path to obtain the secure package. The key which is generated at the time of the job would be released based on your facilities rules for release of PHI. (Virtual media supports Windows®7, 8+, 10)

Epson PP100ii Discproducer

Power: 100VAC - 240VAC ±10%, 50/60Hz
Operating Temp: Temperature: 60° - 86°F (16°C - 30°C)
Physical Specifications: 14.85" W x 15.12" H x 19.41" D
(377mm W x 384mm H x 493mm D) 52.92 Lbs. (24 Kg)

DMC-EZ PC

Power: 100 - 240VAC ±10%, 50/60Hz
Operating Temp: Temperature: 60° - 86°F (16° - 30°C)
Network: 100/1000 Ethernet interface
Physical Specifications:
6.1" W x 13.8" H x 10.8" D (154mm W x 350mm H x 274mm D)
17.62Lbs. (8 Kg)