

- **New RDM tool for interactive, fast and customizable statistical analysis.**
- **Patient dose data (visualization of possible radiation protection problems coming from a modality, drifts in procedures, etc.).**
- **Data related to the facility's service (number of patients, distribution of requests, waiting time between each patient, etc.).**

The dashboard provides a comprehensive overview of radiation dose monitoring. Key features include:

- Population Analysis:** Visualizes patient age distribution and provides summary statistics for total patients, males, and females.
- Modality Overview:** Tracks performance across various modalities, including patient counts, dose rates, and compliance metrics.
- Procedures - Details:** Offers granular data on individual procedures, including patient history, procedure description, and associated radiation doses.
- Performance Overview:** Utilizes heatmaps and detailed tables to analyze the efficiency and radiation output of different procedure types.