




One Platform. Every Format

DMC-EZ 3600

DMC-EZ 3600 DICOM Publishing system, powered by the DMC-EZ Vital Images platform, is a browser-based solution designed to simplify DICOM sharing across your enterprise. Its intuitive **Search-Select-Share** workflow streamlines image access, routing, and delivery while maintaining full compliance with healthcare security and interoperability standards.



Supports: CD/DVD/Blu-ray
 Drives: Two
 Input Bins: Three
 Speed Capability: 30 Full CDs per hour

Robot Options:
 Catalyst II 6000 - 4002359
 Catalyst II 6000N- 4002349

DMC-EZ Vital Images Platform Features:

User Interface & Workflow:

- Intuitive browser-based interface—no local software installation required.
- EZDICOM Sender enables import of external images from any location on the hospital network.

Image Sharing Options

- CD/DVD Creation – Standard media sharing with labeling support.
- Virtual Media™ Cloud sharing – AWS GovCloud Hosting.
- USB Creation – Automated labeling (add-on hardware/software module).

Automation & Scheduling

- Automatic CD/DVD and USB creation via DICOM Push from PACS.
- Rule-based scheduling with multi-AE Title support for flexible automation.
- Automatic or manual encryption using HIPAA/HITECH-compliant 256-bit AES / FIPS 140-2 standards.

Reports & Viewers

- Include reports via EPIC, Sectra, AGFA, GE RIS and Centricity, McKesson, Intelispace, HL7, SR, or Mitra Broker.
- Windows/macOS viewers and built-in web viewer.

Standards & Integration

- IHE DICOM-compliant CD/DVD, USB, and Cloud.
- REST API integration with HIS/RIS/EMR.

Administration & Compliance

- Full HIPAA audit logging.
- LDAP login integration.
- Multilingual support: EN, FR, IT, PT, ES.
- Patented technology.

309 Laurelwood Rd. Suite 20
 Santa Clara, CA 95054
 sales@sstgroup-inc.com

Sales Support: 800.944.6281
 Outside US: 408.350.3450
 Fax: 408.350.3100