



309 Laurelwood Rd. Suite 20
 Santa Clara, CA 95054
 sales@sstgroup-inc.com

Sales Support: 800.944.6281
 Outside US: 408.350.3450
 Fax: 408.350.3100

DMC - EZ DICOM Printer-LF

DICOM Image Distribution

The DMC-EZ DICOM printer solution is a cost-effective way for your institution to provide reference hard copies of studies for your patients and referring physicians and makes it easy to output images directly to paper, replacing expensive film.

DMC-EZ DICOM Printer software functions as a DICOM print server and enables you to send images/studies via DICOM Print. The result is near photographic-quality color or gray scale images on plain paper.



DMC-EZ DICOM Printer Software

- Supports multiple associations
- Multiple DICOM Print AE Titles (Alias)
- Supports multiple laser and/or inkjet printers
- Customizable header and footers
- DIMSE services to act as SCP
- Configurable support for multiple film to paper sizes (3.5" x 5", 4" x 6", 5" x 7", 8" x 10", 8.5" x 11", 8.5" x 14", A3, A4, A5, A6, B4, B5, half letter, executive, 11" x 17", 13" x 19")
- Maximum print resolution, 4800 x 1200 DPI.
- Expandable with 3 additional 500 sheet paper trays
- Patented technology.

Epson WP-8090

The WorkForce Pro WF-8090 A3 printer, powered by PrecisionCore®, offers 13" x 19" printing. It features color and black print speeds of 24 ISO ppm, a high monthly duty cycle, and saves up to 40 percent on printing costs vs. color lasers. Ultra high-yield cartridges print up to 7000 color pages, with a total paper capacity up to 1830 sheets. And, PCL/genuine Adobe® PostScript® 3™ support.

DMC-EZ DICOM Printer Server and Software

High performance CPU contains the DICOM print server software. DICOM Printer Server PC running Windows® 10 controls WP-8090 and connects to DICOM network.

Epson WF-8090

Wide format printer up to 13" x 19" (A3) using 4-color DURABrite® Pro pigment ink. With ISO print speeds of 24 ISO ppm (black) and 24 ISO ppm (color)

Power: 100VAC - 240VAC ±10%, 50/60Hz
 Operating Temp: Temperature: 50° - 95°F (10°C - 35°C)
 Physical Specifications: 22.3" W x 22.1" H x 35.5" D (566.5mm W x 561.4mm H x 901.7mm D)
 57.1 Lbs. (25.9 Kg)

DMC-EZ DICOM Printer Server

Power: 100 - 240VAC ±10%, 50/60Hz
 Operating Temp: Temperature: 60° - 86°F (16° - 30°C)
 Network: 100/1000 Ethernet interface
 Physical Specifications:
 5.16" W x 2.07" H x 5.16" D (131mm W x 52.5mm H x 131mm D)
 1.3Lbs. (.59 Kg)

