



309 Laurelwood Rd. Suite 20
Santa Clara, CA 95054
sales@sstgroup-inc.com

Sales Support: 800.944.6281
Outside US: 408.350.3450
Fax: 408.350.3100

DMC-EZ USB

The DMC-EZ USB DICOM Publishing system is perfect choice for Large Hospitals and Imaging Centers. The DMC-EZ USB allows users to create DICOM USB for distribution to patients, referring physicians and centers of excellence. The device addresses the need for patients, physicians and facilities who no longer have access to a CD or DVD drive on their workstation or PC. USB is also a convenient media for patients with acute or chronic conditions needing to transport large image data sets to and from specialty providers. The DMC-EZ systems satisfy the requirements of all radiology and PACs systems while giving users an intuitive system that is robust and easy to use.

DMC-EZ Medical Software

- Intuitive browser-based user interface.
- Compatible with large format studies (TOMO, Cardio and Ultrasound)
- Creates IHE DICOM-compliant USB drives
- On-demand write-once USB creation for DICOM studies (Query/Retrieve).
- On-demand cloud delivery of DICOM studies (HITECH compliant) - optional module.
- Automatic write-once USB creation using powerful rule-based scheduling utilizing multiple AE Titles.
- Includes multiple viewers for distribution. Windows viewer supports XP through Windows 10.
- Includes diagnostic reports on USB via HL7 feed, SR, Mitra Broker. We also integrate with major RIS/PACS OEM to obtain diagnostics reports. Contact sales for the options.
- View complete report on USB with images and hyper-links.
- Flexible USB labeling options.
- Study anonymization.
- View thumbnail images prior to disc creation.
- Built-in web server and web viewer for alternative distribution of images and reports.
- Optional REST API allows burner to be integrated into your workflow.
- Automated or Manual USB encryption (HIPAA/HITECH 256 AES/FIPS)
- Manual and Automated Rules-based DICOM routing.
- Robot supplies status monitor
- All activity recorded in audit logs (HIPAA Compliance)
- Back-up configuration and HIPAA logs to USB or network location to allow easy restore.
- Import data from CDs and DVDs your facility receives (EzDicom Sender). - (additional CD/DVD reader required)
- Multiple language support - interface and viewer (English, French, Italian, Portuguese, Spanish).
- User login supports active directory.
- **Patented technology.**

DICOM Image Distribution



Rimage Maestro USB Publishing System

The Maestro system is the first fully automated USB system for printing and recording. Engineered for reliable robust performance with industrial grade robotics, and is designed for high-volume and uninterrupted USB publishing. Integrated DiscWatch status light, Multi line LCD status display. Includes display, keyboard and mouse. Automatically burns and prints up to 100 USB unattended. (Depending on configuration.)

DMC-EZ DICOM Server and Software

High performance PC contains the burning and printing software, DMC-EZ application and image cache. DMC-EZ PC running Windows® 10 Pro controls Maestro and connects to DICOM network or DICOM workstation. Includes display, keyboard and mouse.

Color Thermal Printing on USB Media

The state-of-the-art Rimage Everest™ Everest EncoreCP thermal retransfer printer permanently fuses waterproof, scratch-proof, photo-realistic images to the surface of the USB.

Fully Automated Disc Creation

Includes four 25 USB input magazines and 200 USB output bin.

USB capacities:
8GB, 16GB, 32GB

Cases:
Thumb, Pocket

