




Your DICOM distribution experts.

## DVM-EZ Virtual Media

Introducing DVM-EZ Virtual Media (Zero Trust Image Sharing) a revolution in patient care: Zero Trust Image Sharing redefines the way medical images and results are accessed and shared. With our cutting-edge technology, sending images electronically from a single user interface has never been easier.

Our solution is not only efficient but also incredibly user-friendly. EZ to implement and EZ to use, Zero Trust Image Sharing sets a new standard in patient-centric healthcare. Join the future of medical imaging and experience the ease and convenience firsthand.



### DVM-EZ Medical Software Features:

- Intuitive browser-based user interface. **Simple workflow: Search, Select and Distribute!**
- Import outside images from anywhere on the hospital network using the EZDicom Sender.
- On-demand distribution of DICOM studies (Query/Retrieve) from any PACS/VNA using:  
Cloud distribution- Virtual Media (Zero Trust Image Sharing)
- Automatic or Manual encryption (HIPAA/HITECH 256-bit AES/FIPS 140-2) for Virtual Media.
- Includes multiple viewers for distribution (MAC OS 10.15.7+/Windows viewer supports Windows 7 through Windows 11).
- Automatic or Manual anonymization for Virtual Media.
- Automatic or Manual inclusion of diagnostic reports on using HL7 feed, SR, Mitra Broker optional fetcher for EMR/RIS/PACS OEM
- View complete report on Virtual Media with images and hyper-links with our proprietary windows viewer.
- Compatible with large format studies (TOMO, Cardio and Ultrasound).
- IHE DICOM compliant Virtual Media.
- Built-in web viewer for reference internal visualization of images and reports.
- Manual and Automated Rules-based DICOM routing.
- REST API allows burner to be integrated into your HIS/RIS/EMR.
- Burner supplies and status monitor.
- All activity recorded in audit logs (HIPAA Compliance)
- Back-up configuration and HIPAA logs to network location to allow easy restore.
- Multiple language support - interface and viewer (English, French, Italian, Portuguese, Spanish).
- User login supports active directory integration. LDAP/LDAPS
- Hosting Site Upgrade to AWS GovCloud - Additional Fee
- **Patented technology.**

309 Laurelwood Rd. Suite 20  
Santa Clara, CA 95054  
sales@sstgroup-inc.com

Sales Support: 800.944.6281  
Outside US: 408.350.3450  
Fax: 408.350.3100