



Your DICOM distribution experts.

309 Laurelwood Rd. Suite 20
Santa Clara, CA 95054
sales@sstgroup-inc.com

Sales Support: 800.944.6281
Outside US: 408.350.3450
Fax: 408.350.3100

www.sstgroup-inc.com

DICOM Image Distribution

DVM-EZ DICOM Virtual Media

Zero Trust Sharing of patient images and results

Send images electronically from a single device and user interface.

Patients no longer have to wait at the hospital or emergency department for a disc.

No return trips to the hospital are necessary for patients to pick up a disc.

Images and results are simply emailed to patients from an easy to use interface.

Patients can email their images and results to referring physicians or create their own discs at home.

Stop spending money on FedEx or courier services.

Reduce the need to handle discs and sharing between staff, patients and referring physicians.

EZ to implement and EZ to use.

The DVM-EZ™ DICOM virtual media creator provides your facility with a simple, fast and secure way to get DICOM images and results to your patients, referring physicians and other facilities.

DVM-EZ Medical Software

- Intuitive browser-based user interface.
- Compatible with large format studies (TOMO, Cardio and Ultrasound)
- On-demand cloud delivery of DICOM studies (HITECH compliant).
(Zero Trust Image Sharing)
- Includes Windows viewer that supports Windows® 7 through Windows® 11.
- Include diagnostic reports via HL7 feed, SR, Mitra Broker. We also integrate with major RIS/PACS OEM to obtain diagnostics reports. Contact sales for the options.
- View complete report with images and hyper-links.
- Flexible email template.
- Study anonymization.
- Built-in web server and web viewer for alternative distribution of images and reports.
- Virtual media encryption (HIPAA/HITECH 256 AES/FIPS)
- All activity recorded in audit logs (HIPAA Compliance)
- Back-up configuration and HIPAA logs to USB, CD/DVD or network location to allow easy restore.
- Import data from CDs and DVDs your facility receives (EzDicom Sender).
- Multiple language support - interface and viewer (English, French, Italian, Portuguese, Spanish).
- User login supports active directory.
- Can be virtualized in your IT department
- Patented technology.



DVM-EZ Server and Software

High performance CPU running Windows® 10 Pro or 11 Pro contains the DVM-EZ application and image cache. Connects to DICOM network and/or DICOM workstation(s).

Zero Trust Image Sharing

Simple, fast and secure method to get DICOM images to your patients, referring physicians and other facilities. There are **no per study fee's, no special sign-up or accounts to set-up.** By selecting the study(s), entering an email address and pressing "OK" your facility can send and store up to 400 (1 GB) studies in the cloud at any given time. Each with it's own expiration date. No management of the cloud based data is required. Recipient receives an email with a download path to obtain the secure package. The key which is generated at the time of the job would be released based on your facilities rules for release of PHI.

(Virtual media supports Windows® 7, 8, 10, 11)

DVM-EZ PC

Power: 100 - 240VAC ±10%, 50/60Hz
Operating Temp: Temperature: 60°- 86°F
Network: 100/1000 Ethernet interface
Physical Specifications: 6.1" W x 12.8" H x 11.5" D, 14.59 lbs.